


Convert Your Microscope into a Digital Scanner

MEIJI TECHNO



Efficient manual scanning for any microscope

DSS Basic is a manual whole slide imaging solution that allows users to digitize slides using an existing trinocular microscope.

Operator manually moves the sample, while the software automatically captures and stitches the image into a large digital slide.

Install DSS Basic
on your microscope!

Kit includes:

1. Digital camera (3 MP, 120 FPS)
2. DSS Basic measurement software
3. Your microscope




Quick, easy, and complete digital slide capture


How it works



1 Install the camera




2 Open the DSS Basic software




3 Place the slide on the stage



4 Press "Start" and move the sample




5 Software automatically captures and stitches the field




6 Save and view the digital slide


Digital slide analysis



Database for archive management
Organize, store, and quickly access your scanned images and data.



Count and classification
Automatic or manual object identification on digital slides. Measurement of particle linear and optical parameters.



Marks and comments
Make your comments and put marks on a digital slide.

MEIJITECHNO



MEIJI TECHNO CO., LTD.

322-1, Chikumazawa Miyoshi machi,
Iruma-gun Saitama 354-0043, Japan
tel.: (0)49-259-0111
fax: (0)49-259-0113
meiji@meijitechno.co.jp
www.meijitechno.co.jp



MEIJI TECHNO AMERICA

591 Division St., Campbell CA 95008
Toll free USA: 800-832-0060
tel.: 408-226-3454
fax: 408-226-0900
admin@meijitechno.com
www.meijitechno.com



MEIJI TECHNO EUROPE

Brown-Boveri-Straße 6, B17-1
2351 Wiener Neudorf, Austria
tel.: +43 (0) 2236 892465
fax: +43 (0) 2236 892464
info@meijitechno-europe.com
www.meijitechno-europe.com



MEIJI TECHNO UK, LTD.

Meiji Techno UK Limited,
The Vineyard, Axbridge
BS26 2AW, SOMERSET
tel.: 01934 733655
enquiries@meijitechno.co.uk
www.meijitechno.co.uk

| www.meijitechno.com |

We reserve the right to change specification without notice.

Official Distributor

Rev 1.0_03.2026 EN